

# Stock Solution



TD-S Revision 2.0

Creation Date:  
Revision Date: 10/17/2019

## 25-50 mg/ml Chloramphenicol Stock Solution

### Instructions

1. Weigh 0.25-0.5 g of Chloramphenicol ([Chloramphenicol, GoldBio Catalog # C-105](#) [CAS 56-75-7, mw. = 323.1]).
2. Add 10 ml of 95% EtOH (70% may also be used). Dissolve completely.
3. Stock may be kept at -20°C for 1 year.