

Stock Solution



TD-S Revision 2.0

Creation Date: 5/24/2017
Revision Date: 10/21/2019

100 mg/ml Vancomycin Stock Solution

Instructions

1. Weigh 1 g of Vancomycin ([Vancomycin Hydrochloride, GoldBio Catalog # V-200](#), [CAS 1404-93-9, mw. = 1485.71 g/mol]).
2. Add 10 ml of sterile H₂O. Dissolve completely.
3. Prewet a 0.22 µm syringe filter by drawing through 5-10 ml of sterile H₂O and discard water.
4. Sterilize Vancomycin Stock through the prepared 0.22 µm syringe filter.
5. Stock may be kept at -20°C for 1 year.