Stock Solution



TD-S Revision 2.0

Creation Date: 10/28/2015 Revision Date: 10/22/2019

50 mg/ml Aztreonam, Solubilized Stock Solution

Instructions

- 1. Weigh 0.50 g of Solubilized Aztreonam (<u>Aztreonam, Solubilized, GoldBio Catalog # A-940</u> [CAS 78110-38-0, mw. = 435.43 g/mol]).
- 2. Add 10 ml of sterile H₂O. Dissolve completely.
- 3. Prewet a 0.22 μm syringe filter by drawing through 5-10 ml of sterile H_2O and discard water.
- 4. Sterilize Aztreonam stock through the prepared 0.22 μm syringe filter.
- 5. The stability of Aztreonam in solution is not known so we recommend preparing it fresh. Stock solutions of some antibiotics are often kept in small aliquots at -70°C for as long as 3-12 months.