

Stock Solution



TD-S Revision 2.0

Creation Date: 6/26/2014
Revision Date: 10/17/2019

10mg/ml Bialaphos Stock Solution

Instructions

1. Weigh 0.1 g of Bialaphos ([Bialaphos, GoldBio Catalog # B0178](#) [CAS 71048-99-2, mw. = 345.26]).
2. Add 10 ml of sterile H₂O. Dissolve completely.
3. Prewet a 0.22 µm syringe filter by drawing through 5-10 ml of sterile H₂O and discard water.
4. Sterilize Bialaphos Stock through the prepared 0.22 µm syringe filter.
5. Stock may be kept at -20°C for 1 year.