

Stock Solution



TD-S Revision 2.0

Creation Date: 7/18/2014

Revision Date: 4/16/2019

L-Arabinose 10% w/w Stock Solution

Instructions

1. Add 1g of [L-Arabinose \(GoldBio Catalog # A-300\)](#) to 9ml of milli-Q water in a small falcon tube.
2. Add milli-Q water to bring volume up to 10 ml.
3. Filter sterilize with a 0.22m filter.
4. Can be stored at room temperature.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com