

# Stock Solution



TD-S Revision 2.0

Creation Date: 7/18/2014

Revision Date: 4/12/2019

## Kinetin 1mM Stock Solution

### Instructions

1. Dissolve 21.5 mg of [Kinetin \(GoldBio Catalog # K-100\)](#) in 1 ml of 1N NaOH.
2. Bring solution to 10 ml with molecular biology grade water to get a final solution of 10mM.
3. Sterilize by filtration.
4. Store Kinetin Stock Solution at -20°C.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: [www.goldbio.com](http://www.goldbio.com)

Email: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)