

Stock Solution



TD-S Revision 2.0

Creation Date:
Revision Date: 10/22/2019

1-5 mg/ml Clindamycin Stock Solution

Instructions

1. Weigh 0.01-0.05 g of Clindamycin ([Clindamycin Hydrochloride, GoldBio Catalog # C-175](#), [CAS 21462-39-5, mw. = 479.5 g/mol]).
2. Add 10 ml of sterile H₂O. Dissolve completely
3. Prewet a 0.22 μm syringe filter by drawing through 5-10 ml of H₂O and discard water.
4. Sterilize Clindamycin Stock through the prepared 0.22 μm syringe filter.
5. Stock may be kept at -20°C for 1 year.