

# Stock Solution



TD-S Revision 2.0

Creation Date: 7/21/2015

Revision Date: 4/12/2019

## 10% CHAPS – 50 ml Stock Solution

### Instructions

1. Weigh 5 g of CHAPS ([CHAPS, Gold Biotechnology, Catalog# C-080](#) [CAS 75621-03-3, mw. = 614.88 g/mol]).
2. Fill to final volume of 50 ml with sterile dH<sub>2</sub>O.
3. Store CHAPS Stock Solution at room temperature.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: [www.goldbio.com](http://www.goldbio.com)

Email: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)