

Stock Solution



TD-S Revision 2.0

Creation Date: 10/15/2015

Revision Date: 10/22/2019

100mM AEBSF Stock Solution

Instructions

1. Weigh 0.9592 g of AEBSF ([AEBSF, GoldBio Catalog # A-540](#) [CAS 30827-99-7, mw. = 239.7]).
2. Fill to a final volume of 40 ml with pure sterile H₂O.
3. Mix until completely dissolved.
4. Use immediately or store in small aliquots at -20°C for months.

Note: AEBSF is effective at a concentration of 0.1-1mM. Dilute the stock solution 1:100 for a 1mM working concentration. AEBSF is also soluble ethanol.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com