

Stock Solution



TD-S Revision 2.0

Creation Date: 10/15/2015

Revision Date: 10/17/2019

10 mg/ml Aprotinin Stock Solution

Instructions

1. Weigh 25 mg of Aprotinin ([Aprotinin, GoldBio Catalog # A-655](#) [CAS 9087-70-1, mw. = 6511.44 g/mol]).
2. Fill to a final volume of 2.5 ml with Phosphate Buffered Saline ([PBS \(Phosphate Buffered Saline\) Tablets, GoldBio Catalog # P-271](#)).
3. Mix until completely dissolved.
4. Use immediately or store in small aliquots at -20°C for months.

Note: Aprotinin is effective at a concentration of 0.6-4 µg/ml. Dilute the stock solution 1:5000 for a working concentration of 2 µg/ml.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com