

Stock Solution



TD-S Revision 2.0

Creation Date: 7/10/2015

Revision Date: 4/12/2019

50 mg/ml Azlocillin Sodium Stock Solution

Instructions

1. Weigh 0.50 g of Azlocillin Sodium Salt ([Azlocillin Sodium Salt, GoldBio Cat# A-920](#) [CAS 37091-65-9, mw. = 483.47 g/mol]).
2. Add 10 ml of sterile H₂O. Dissolve completely.
3. Prewet a 0.22 µm syringe filter by drawing through 5-10 ml of sterile H₂O and discard water.
4. Sterilize Azlocillin stock through the prepared 0.22 µm syringe filter.
5. Stock may be stored in small aliquots at -70°C for up to 1 year.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com