

# Stock Solution



TD-S Revision 2.0

Creation Date: 8/29/2014  
Revision Date: 10/17/2019

## 50 mg/ml Erythromycin Stock Solution

### Instructions

1. Weigh 0.5 g of Erythromycin [[Erythromycin, GoldBio Catalog # E-350](#) (CAS 114-07-8, MW = 733.93)].
2. Add 10 ml of 95% EtOH (lower concentrations may be used). Dissolve completely.
3. Stock may be kept at -20°C for 1 year.