

# Stock Solution



TD-S Revision 2.0

Creation Date: 4/14/2015

Revision Date: 4/19/2019

## 1% Ponceau S Stock Solution

### Instructions

1. Weigh 1 g of Ponceau S ([Ponceau S, GoldBio Catalog # P-330](#) [CAS 6226-79-5, mw. = 760.57]).
2. Add 50 ml Acetic Acid. Dissolve completely.
3. Fill to 1 L with distilled H<sub>2</sub>O.
4. Stock solution may be stored at room temperature or 4°C. Do not freeze.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: [www.goldbio.com](http://www.goldbio.com)

Email: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)