

Stock Solution



TD-S Revision 2.0

Creation Date: 10/6/2015
Revision Date: 10/23/2019

Nisin Stock Solution – 1 mg/ml

Instructions

1. Weigh 1 mg of Nisin ([Nisin, GoldBio Catalog # N-660](#) [CAS 1414-45-5, mw. = 3,354.07 g/mol]) into a 1.5 ml centrifuge tube.
2. Add 1 ml of 0.05% acetic acid.
3. Mix for 5-10 minutes at room temperature.

Note: In Nisin, insoluble proteins will form a cloudy suspension.