

# Stock Solution



TD-S Revision 1.0

Creation Date: 8/6/2021

Revision Date: 8/6/2021

## 100 mg/ml Abscisic Acid Stock Solution

### Instructions

1. Weigh 1 g of Abscisic Acid ([Abscisic Acid, GoldBio Catalog # A-050](#) [CAS 21293-29-8, mw. = 264.32]).
2. Add 10 ml of 95% EtOH (DMF or DMSO may also be used). Dissolve completely.
3. Stock may be kept at -20°C for 1 year.

Alternatively, Abscisic Acid is soluble in 0.1N NaOH at 20 mg/ml.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: [www.goldbio.com](http://www.goldbio.com)

Email: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)