

Stock Solution



TD-S Revision 1.0

Creation Date: 12/1/2021

Revision Date: 12/1/2021

Lysozyme Stock Solution

Instructions

1. Weigh out 10 mg of Lysozyme, Egg White into a microcentrifuge tube.
2. Add 1 mL of 10mM Tris-HCl, pH 8.0 and mix by pipet.
3. Store at 4°C for up to 1 month.

Gold Biotechnology

St. Louis, MO

Ph: (314)890-8778

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com