

# Stock Solution



TD-S Revision 1.0

Creation Date: 5/27/2022

Revision Date: 6/8/2022

## 1 mg/ml (2.55mM) Dexamethasone Stock Solution

### Instructions

1. Weigh 5 mg of dexamethasone ([Dexamethasone, GoldBio Catalog # D-940](#)).
2. Add 5 mL ethanol.
3. Vortex until complete solubilization.
4. Aliquot and store at 4°C.

Once dexamethasone is solubilized, dilutions can be prepared by adding sterile culture medium.

Gold Biotechnology

St. Louis, MO

Ph: (800) 248-7609

Web: [www.goldbio.com](http://www.goldbio.com)

Email: [contactgoldbio86@goldbio.com](mailto:contactgoldbio86@goldbio.com)