

**CROSSLINKED AGAROSE BEADS
 TECHNICAL SPECIFICATIONS
 Crosslinked Agarose Bead Resins**

BEAD (Geometry)	Spherical				
CROSS-LINKED	Yes				
CHEMICAL STABILITY	Stable in very strong acid and basic solution and also in dissociation reagents				
MECHANICAL STABILITY	Suitable and comparable to commercial counterparts				
AGAROSE %	2%	4%		6%	
BEAD SIZE RANGE	Standard	Standard	Fine	Standard	Fine
	50-150 µm	50-150 µm	20-50 µm	50-150 µm	20-50 µm
CATALOG #	A-191	A-192	A-171	A-193	A-172
ANTIMICROBIAL AGENT	20% Ethanol				
STORAGE TEMPERATURE	2 - 25°C				

For laboratory use only. Not for use in diagnostic or therapeutic procedures.