

Stock Solution



TD-S Revision 2.0

Creation Date: 8/12/14
Revision Date: 10/21/2019

Preparation of Salmon Gal

Instructions

1. To prepare a stock solution of Salmon Gal ([Salmon Gal, GoldBio Catalog # C-255](#), [CAS 138182-21-5, mw. = 371.4 g/mol]) add 25 mg of Salmon Gal to 1 mL of dimethylformamide (DMF) or dimethylsulfoxide (DMSO) for a final concentration of 25 mg/ml.
2. There is no need to sterilize the stock solutions when the solution is made with DMF or DMSO.
3. The recommended working concentration for Salmon Gal is 150 µg/ml.