

Growth Factor Data Sheet

GoldBio growth factors are manufactured for **RESEARCH USE ONLY** and cannot be sold for human consumption!

IL1B is a pro-inflammatory cytokine produced in a variety of cells including monocytes, tissue macrophages, keratinocytes and other epithelial cells. Both IL1A and IL1B bind to the same receptor and have similar, if not identical, biological properties. These cytokines have a broad range of activities including, stimulation of thymocyte proliferation by inducing IL2 release, B-cell maturation and proliferation, mitogenic FGF-like activity and the ability to stimulate the release of prostaglandin and collagenase from synovial cells. However, whereas IL1A is predominantly a cell associated cytokine, IL1B is a secreted cytokine. The 17.3 kDa mature rhesus IL1B shares 98% amino acid sequence identity with human IL-1B.

Catalog Number	1210-01B
Product Name	IL1B (IL-1 beta), Rhesus macaque Recombinant Rhesus macaque Interleukin-1 β IL1 β IL 1-Beta Interleukin1B
Source	<i>Escherichia coli</i>
MW	~17.3 kDa (153 amino acids)
Sequence	APVRS LHCTL RDAQ LKSLVM SGPYELKALH LQGQDLEQQV VFSMSFVQGE ESNDKIPVAL GLKAKNLYLS CVLKDDKPTL QLESVDPKNY P KKKMEKRFV FNKIEINNKL EFESAQFPNW YISTSQAENM PVFLGGTRGG QDITDFTMQF VSS
Accession Number	P48090
Purity	>98% by SDS-PAGE and HPLC analyses
Biological Activity	Fully biologically active when compared to standard. The ED ₅₀ as determined by a cell proliferation assay using murine D10.G4.1 cells is less than 10 pg/ml, corresponding to a specific activity of >1.0 × 10 ⁸ IU/mg.
Formulation	Sterile filtered white lyophilized powder. Purified and tested for use in cell culture.
Storage/Handling	This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage. The reconstituted sample can be apportioned into working aliquots and stored at -80 °C for up to 6 months. Avoid repeated freeze/thaw cycles.
Reconstitution	The sample should be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in a siliconized tube using PBS that contains a 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Reconstituted solutions are stable for up to one week at 2-8°C. Stock solutions should be aliquoted and stored at -80°C. Further dilutions should be made in appropriate buffered solutions containing BSA or serum.