

Stock Solution



TD-S Revision 2.0

Creation Date: 8/1/2014
Revision Date: 10/22/2019

Luminol (Sodium Salt) Stock Solution

Instructions

1. Dissolve luminol sodium salt ([GoldBio Catalog # L-145](#), [CAS 20666-12-0, mw. = 199.150 g/mol]) powder in PBS to a final concentration of 50 mg/ml.
2. This solution should be filtered through a 0.2 μm filter.
3. Aliquot into sterile Eppendorf tubes.
4. Store at -20°C .