

**PLAIN AGAROSE BEADS  
 TECHNICAL SPECIFICATIONS  
 Plain Agarose Bead Resins**

<b>BEAD (Geometry)</b>	Spherical						
<b>CROSS-LINKED</b>	No						
<b>CHEMICAL STABILITY</b>	Stable in strong acid and basic solutions						
<b>MECHANICAL STABILITY</b>	Suitable and comparable to commercial counterparts						
<b>AGAROSE %</b>	2%	4%		6%		8%	10%
<b>BEAD SIZE RANGE</b>	Standard	Standard	Fine	Standard	Fine	Standard	Standard
	50-150 µm	50-150 µm	20-50 µm	50-150 µm	20-50 µm	50-150 µm	50-150 µm
<b>CATALOG #</b>	<a href="#">A-181</a>	<a href="#">A-182</a>	<a href="#">A-173</a>	<a href="#">A-183</a>	<a href="#">A-174</a>	<a href="#">A-184</a>	<a href="#">A-185</a>
<b>ANTIMICROBIAL AGENT</b>	20% Ethanol						
<b>STORAGE TEMPERATURE</b>	2 - 8°C						

**For laboratory use only. Not for use in diagnostic or therapeutic procedures.**