

# Growth Factor Data Sheet

GoldBio growth factors are manufactured for **RESEARCH USE ONLY** and cannot be sold for human consumption!

Murine CCL8 is a chemokine containing four conserved cysteine residues, the first two of which are adjacent. It is expressed by peripheral blood leukocytes, as well as a variety of other cells. Expression of CCL8 is induced by IFNG, IFNB, IL1B, and various cytokine inducers such as viruses, endotoxins, and mitogens. It is chemotactic for monocytes, lymphocytes, basophils, and eosinophils and plays a role in tumor-associated leukocyte infiltration by recruiting leukocytes to sites of inflammation. CCL8 also inhibits CCR5-mediated HIV-1 entry and replication. CCL8 is a ligand for CCR1, CCR2B, CCR3, CCR5, CCR8, ACKR2, and ACKR4.

<b>Catalog Number</b>	<b>2340-08</b>
<b>Product Name</b>	<b>CCL8, Murine</b> Recombinant Human Chemokine (C-C motif) Ligand 8 (CCL8) Monocyte Chemoattractant Protein 2 (MCP2, MCP-2) Small Inducible Cytokine Subfamily A Member 8 (SCYA8) HC14
<b>Source</b>	<i>Escherichia coli</i>
<b>MW</b>	~8.5 kDa (74 amino acids)
<b>Sequence</b>	<a href="#">Q9Z121</a>
<b>Accession Number</b>	GPDKAPVTCC FHVLLKLIPL RVLKSYERIN NIQCPMEAVV FQTKQGMSLC VDPTQKWWSE YMEILDQKSQ ILQP
<b>Purity</b>	>97% by SDS-PAGE and HPLC analyses
<b>Biological Activity</b>	Fully biologically active when compared to standard. The biological activity determined by a chemotaxis bioassay using human peripheral blood monocytes is in a concentration range of 10-100 ng/ml.
<b>Formulation</b>	Sterile filtered white lyophilized powder. Purified and tested for use in cell culture.
<b>Storage/Handling</b>	This lyophilized preparation is stable at 2-8°C, but should be kept at -20°C for long term storage. The reconstituted sample can be apportioned into working aliquots and stored at -80 °C for up to 6 months. Avoid repeated freeze/thaw cycles.
<b>Reconstitution</b>	The sample should be briefly centrifuged prior to opening to bring the contents to the bottom. Reconstitute in a siliconized tube using PBS that contains a 0.1% BSA to a concentration of 0.1-1.0 mg/mL. Reconstituted solutions are stable for up to one week at 2-8°C. Stock solutions should be aliquoted and stored at -80°C. Further dilutions should be made in appropriate buffered solutions containing BSA or serum.