

**AFFINITY HIS-TAG PURIFICATION  
TECHNICAL SPECIFICATIONS  
Nickel NTA Agarose BULK resins**

<b>BEAD (Geometry, size)</b>	Spherical, Standard: ~ 50-170 µm
<b>CROSS-LINKED</b>	Yes
<b>AGAROSE %</b>	6% agarose
<b>LIGAND</b>	Nitrilotriacetic Acid (NTA)
<b>PRODUCT NAME</b>	Nickel NTA Agarose Resin
<b>Catalog Number</b>	<a href="#">H-350</a>
<b>BINDING/LOADING CAPACITY</b>	>50 mg/ml gel
<b>ANTIMICROBIAL AGENT</b>	15% Ethanol
<b>STORAGE TEMPERATURE</b>	2 - 8°C

**For laboratory use only. Not for use in diagnostic or therapeutic procedures.**