

Stock Solution



TD-S Revision 2.0

Creation Date: 7/7/2015
Revision Date: 10/21/2019

25 mg/ml Trimethoprim Lactate Stock Solution

Instructions

1. Weigh 0.25 g of Trimethoprim Lactate ([Trimethoprim Lactate, GoldBio Catalog # T-351](#) [CAS 23256-42-0, mw. = 380.40 g/mol]).
2. Add 10 ml of DMSO ([Sterile Filtered DMSO, Ultra Pure, GoldBio Catalog # D-361](#) [CAS 67-68-5, mw. = 78.13 g/mol]). Dissolve completely.
3. Stock may be kept at -20°C for up to 1 year.