

**Streptavidin Agarose Resin
TECHNICAL SPECIFICATIONS
Agarose Bead Bulk Resins**

BEAD (Geometry)	Spherical, Standard: ~50-150 μ m
CROSS-LINKED	Highly Crosslinked
FREE BIOTIN BINDING CAPACITY	>120 nmol/ml gel
AGAROSE %	6%
CATALOG #	S-105
ANTIMICROBIAL AGENT	20% Ethanol
STORAGE TEMPERATURE	4°C

For laboratory use only. Not for use in diagnostic or therapeutic procedures.

Gold Biotechnology

St. Louis, MO

Web: www.goldbio.com

Ph: (314) 890-8778

email: contactgoldbio86@goldbio.com