

Stock Solution



TD-S Revision 1.0

Creation Date: 3/12/2019

Revision Date: 3/12/2019

10 mg/ml Bleomycin Sulfate Stock Solution

Instructions

1. Weigh 100 mg of Bleomycin Sulfate ([Bleomycin Sulfate, GoldBio Catalog # B-910](#) [CAS 9041-93-4, mw. = 1,512.62 g/mol]).
2. Add 10 ml of sterile H₂O. Dissolve completely.
3. Prewet a 0.22 µm syringe filter by drawing through 5-10 ml of sterile H₂O and discard water.
4. Sterilize Bleomycin Sulfate Stock through the prepared 0.22 µm syringe filter.
5. Stock may be kept at -20°C for up to 3 months.

Gold Biotechnology

St. Louis, MO

Ph: (314)890-8778

Web: www.goldbio.com

Email: contactgoldbio86@goldbio.com