

Competent Cell Recovery Medium

Introduction

GoldBio's Competent Cell Recovery Medium has a superior formulation recommended for all competent cell transformation. Our Recovery Medium is formulated to achieve the highest possible transformation efficiency from any competent cell.

Materials

- Competent Cell Recovery Medium (GoldBio Catalog # [CC-300](#))

Storage and Handling

- This product may be shipped frozen but should be stored at 4°C immediately upon arrival. When stored under the recommended conditions and handled correctly, these products should be stable for at least 1 year from the date of receipt.

Associated Products

- [DH10B Chemically Competent *E. coli* Cells \(GoldBio Catalog # CC-100\)](#)
- [DH5-alpha Chemically Competent *E. coli* Cells \(GoldBio Catalog # CC-101\)](#)
- [BL21 Chemically Competent *E. coli* Cells \(GoldBio Catalog # CC-102\)](#)
- [BL21 \(DE3\) Chemically Competent *E. coli* Cells \(GoldBio Catalog # CC-103\)](#)
- [DL39 \(DE3\) Chemically Competent *E. coli* Cells \(GoldBio Catalog # CC-104\)](#)
- [DH10B Electrocompetent *E. coli* Cells \(GoldBio Catalog # CC-200\)](#)
- [DH10B-Pro™ Electrocompetent *E. coli* Cells \(GoldBio Catalog # CC-201\)](#)
- [DH10B BAC Library Electrocompetent *E. coli* Cells \(GoldBio Catalog # CC-202\)](#)
- [DH5-alpha Electrocompetent *E. coli* Cells \(GoldBio Catalog # CC-203\)](#)
- [BL21 \(DE3\) Electrocompetent *E. coli* Cells \(GoldBio Catalog # CC-204\)](#)
- [TG1 Phage Display Electrocompetent Cells \(GoldBio Catalog # CC-205\)](#)
- [SS320 Phage Display Electrocompetent Cells \(GoldBio Catalog # CC-206\)](#)
- [GV3101 Agrobacterium Electrocompetent Cells \(GoldBio Catalog # CC-207\)](#)
- [AGL-1 Agrobacterium Electrocompetent Cells \(GoldBio Catalog # CC-208\)](#)
- [Lactococcus lactis MG1363 Competent Cells \(GoldBio Catalog # CC-209\)](#)
- [Lactococcus lactis IL1403 Competent Cells \(GoldBio Catalog # CC-210\)](#)