

Competent Cell Recovery Medium

Introduction

GoldBio's Universal Competent Cell Recovery Medium is specifically formulated for optimal competent cell transformations of every kind of cell, including *E. coli* and *Agrobacterium*.

GoldBio recovery media is included with all of our competent cell kits.

Materials

- Competent Cell Recovery Medium (GoldBio Catalog # CC-300)

Storage and Handling

- This product may be shipped frozen but should be stored at 4°C immediately upon arrival. When stored under the recommended conditions and handled correctly, these products should be stable for at least 1 year from the date of receipt.

Associated Products

- DH10B Chemically Competent *E. coli* Cells (GoldBio Catalog # CC-100)
- DH5-alpha Chemically Competent *E. coli* Cells (GoldBio Catalog # CC-101)
- BL21 (DE3) Chemically Competent *E. coli* Cells (GoldBio Catalog # CC-103)
- DL39 (DE3) Chemically Competent *E. coli* Cells (GoldBio Catalog # CC-104)
- GV3101 *Agrobacterium* Chemically Competent Cells (GoldBio Catalog # CC-105)
- AGL-1 *Agrobacterium* Chemically Competent Cells (GoldBio Catalog # CC-106)
- LBA4404 *Agrobacterium* Chemically Competent Cells (GoldBio Catalog # CC-107)
- EHA105 *Agrobacterium* Chemically Competent Cells (GoldBio Catalog # CC-108)
- HB101 Chemically Competent *E. coli* Cells (GoldBio Catalog # CC-150)
- DH10B Electrocompetent *E. coli* Cells (GoldBio Catalog # CC-200)
- DH10B-Pro™ Electrocompetent *E. coli* Cells (GoldBio Catalog # CC-201)
- DH5-alpha Electrocompetent *E. coli* Cells (GoldBio Catalog # CC-203)
- BL21 (DE3) Electrocompetent *E. coli* Cells (GoldBio Catalog # CC-204)
- GV3101 *Agrobacterium* Electrocompetent Cells (GoldBio Catalog # CC-207)
- AGL-1 *Agrobacterium* Electrocompetent Cells (GoldBio Catalog # CC-208)
- LBA4404 *Agrobacterium* Electrocompetent Cells (GoldBio Catalog # CC-220)
- EHA105 *Agrobacterium* Electrocompetent Cells (GoldBio Catalog # CC-225)